## Dimensions

205mm

# INSTYLAN

INSTRUCTIONS FOR USE

105 mm

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## Sterile solution for intravesical irrigation based on hyaluronic acid

**INSTYLAN****Description**

**INSTYLAN** is a sterile, transparent, colourless viscous solution based on hyaluronic acid of non animal origin with physiological Ph.

**INSTYLAN** is a viscoelastic protective agent of the intercellular substance of the bladder and urethra.

**Classification of the product**

Medical device – sterilized by steam – pyrogen free – Class IIa

**Package**

50 ml pre-filled disposable plastic bag which contains 0.16% (80 mg/ 50 ml) of sodium hyaluronate.

**Each package contains:**

- 1 Polymer bag with a nozzle for the urinary catheter in a carton.

**Quality composition:**

50 ml of sterile solution for intravesical irrigation based on hyaluronic acid **INSTYLAN** 0.16% contains:

Sodium hyaluronate

Excipients: Sodium Chloride ,Water for injection, Phosphate buffer

**Indications**

- Deficiency of a layer of glycosaminoglycans (GAG) and restoration of the protective layer of the bladder mucosa after invasive interventions (such as cystoscopy, transurethral resection (TUR) of the prostate, ureteroscopy, etc.).
- Recurrent bacterial cystitis.
- Interstitial cystitis.
- Painful bladder syndrome.
- Hemorrhagic cystitis.
- Actinic cystitis
- Chemotherapy-induced cystitis.

**Mode of action**

Hyaluronic acid is contained in all tissues and is the most significant component of the intercellular matrix and thus constitutes a significant part of the glycosaminoglycans of the urothelium of the bladder and urethra.

When injected into the bladder cavity, it creates a natural-like, protective visco-elastic layer on the surface of the urothelium, thus protecting the bladder mucosa from irritating and damaging factors.

**Contraindications**

- Do not use the product in patients with hypersensitivity to hyaluronic acid or those with a history of allergic reaction to any of the constituents of solution.
- Do not use the product during pregnancy and lactation.
- Do not use the product in children and adolescents aged under 18 years.

**Method of application**

**INSTYLAN** is intended only for instillation into the bladder. This procedure should be performed by a specialist with the appropriate experience and medical education to comply with all the standards of the urinary catheterization.

**INSTYLAN** is intended for intravesical administration, once weekly. The course is from 4 to 12 instillations.

Before administration of **INSTYLAN**, the bladder must be emptied.

**INSTYLAN** is injected into the bladder cavity by means of an urinary catheter. Carefully follow the operations suggested by the normal protocol for vesical catheter management. After administration of the solution, the catheter is removed.

The recommended retention time of the solution in the bladder lumen is 30 minutes to 2 hours.

Then the bladder should be emptied physiologically.

Prior to use, the bag with the solution should be warmed to body temperature.

The solution warmed to body temperature reduces the patient's discomfort. This procedure should only be carried out using dry heat (for example, heating pad or hot plate).

The solutions should not be heated in water due to increased risk of contamination. To avoid potential damage to the polymer bag and damage or discomfort to the patient, the solutions should not be heated in the microwave.

Strict aseptic rules should be followed during the entire procedure.

**How to use:**

- Verify the integrity of the external packaging, the expiry date, the volume of the device
- Open the plastic bag to take out the plastic device
- Control the integrity of the plastic device
- Read the instruction for use
- Remove the orange cap from the luer lock of the plastic device

- Insert the luer-lock in the urinary catheter.
- Break the internal red cone of the luer-lock to release the solution into urological catheter
- Push the plastic device to allow the exit of entire volume solution.
- Performs the procedure and continuously check visually

#### Adverse reactions

Hypersensitivity reactions. In rare cases: bladder irritation, bladder cramps, dysuria, pain, diarrhea, dizziness. In case of any adverse reactions, you should immediately stop using the medical device and inform the doctor and the manufacturer.

#### Warnings

- The product must not be administrated orally
- Do not use if the solution is not transparent and colorless
- The product must be used immediately after opening
- The product is intended for single use only.
- Do not resterilize. Repeated sterilization may cause cross-contamination
- Do not use if the packaging has been opened or is damaged
- Do not use after the expiry date.



**Shelf life:** 2 years in an intact package.

#### Storage conditions

Store at temperature between +5 ° to 30 °C inclusive, protected from light. Do not freeze. Keep out of reach of children.

#### Disposal

Dispose the materials used according to the applicable laws on medical waste disposal requirements.

 <b>Diaco Biofarmaceutici S.R.L.</b> Via Flavia 124, 34147 Trieste-Italy		Symbols used on the packaging of the box	
			Do not reuse - for single use
	Batch code		Sterilized using steam or dry heat
	Expiry date		Do not resterilize
	Do not use if the package has been opened or is damaged		Protect from light and heat
	Consult instruction for use		Temperature limit for warehouse storage
	Non - pyrogenic		Date of manufacture

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